Gebrauchsanweisung/Instructions for use/Mode d'emploi/
Instrucciones de uso/Istruzioni per l'uso/Instrukcja obsługi

AXION STIM-PRO X9B

Digitaler 4-Kanal Nerven- und Muskelstimulator/Digital 4-channel
nerve and muscle stimulator/Stimulateur numérique à 4 canaux
pour les nerfs et les muscles/Estimulador nervioso y muscular digital
de 4 canales/Elettrostimolatore nervoso e muscolare digitale a
4 canali/4-kanátowy cyfrowy stymulator nerwów i mięśni

INHALTSVERZEICHNIS / CONTENT

2 – 34	Gebrauchsanweisung
35 – 64	Instructions for use
65 – 98	Mode d'emploi
99 – 131	Instrucciones de uso
132 – 166	Istruzioni per l'uso
167 – 199	Instrukcja obsługi

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ALLGEMEINES

Was ist TENS?

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GENERAL INFORMATION

What is TENS?

TENS (transcutaneous electrical nerve stimulation) is a proven procedure in pain therapy. It is used for symptomatic relief of acute, chronic and intractable pain.

Functional principle of TENS

Transcutaneous electrical nerve stimulation is a non-invasive, drug-free method of pain management. TENS uses tiny electrical pulses sent through the skin to the nerves to alter your perception of pain. TENS does not cure a physiological problem; it only helps control pain. TENS doesn't work for everyone, but for most patients it helps reduce or eliminate pain.

What is EMS?

EMS (electrical muscle stimulation) is a method which uses the beneficial effect of low or medium frequency current for the treatment of all muscles.

This method is highly intensive and can support increasing **lean body mass** and muscle building, relieves tension and muscular imbalances.

Functional principle of EMS

In this type of procedure, the muscles are stimulated with electrical impulses. Low or medium frequencies are used, which are safe for humans when used correctly. The electrical impulses trigger a short muscle twitch. Several of these impulses in succession cause a complete muscle contraction.

IMPORTANT NOTICES

- ▶ Read the instructions for use before use. Make sure that you observe all warnings and safety instructions. Ignoring the instructions may result in user injury or device damage.
- ▶ Report all serious incidents that have occurred in connection with the product to the manufacturer and the competent authority of the Member State (where the user is located).
- ▶ STIM-PRO X9B supplies energy to the human body by means of electrodes or—in the case of treatment of incontinence—probe electrodes. CE 2460 only applies to the TENS/EMS device and the probe electrodes. CE is valid for the electrodes (class I products).
- ▶ **Consult a medical professional before use!**

EN

INTENDED USE

The intended purpose of the STIM-PRO X9B is the transcutaneous (through the skin) and transmucosal (through the mucous membrane) electrical nerve and muscle stimulation in humans to relieve acute, chronic and intractable pain, mitigate incontinence, and strengthen muscles. The device has been developed for home and professional use. For home use, application is permitted from the age of 18 only.

INDICATIONS/CONTRAINDICATIONS

Indications

TENS is used to treat the following pain patterns: muscle, joint or bone pain. Pain associated with diseases such as arthrosis or fibromyalgia, conditions such as back pain, neck pain or tendinitis.

EMS is used to relax muscle spasms, prevent or delay muscle atrophy, improve local circulation, restore muscle tone, maintain or increase range of motion, and treat incontinence.

Contraindications

- ▶ Patients with an implanted electronic device, such as a pacemaker or defibrillator, or any other metallic or electronic implant may not undergo treatment without first consulting a physician.
- ▶ Patients with heart disease, epilepsy, cancer, fever or other illnesses may not undergo treatment without prior medical consultation.
- ▶ Use the electrodes only on healthy, uninjured skin.
- ▶ If heart disease is suspected, the treatment should only be performed under medical supervision.
- ▶ Do not stimulate via the carotid sinus nerve, especially in patients with a known sensitivity to the carotid sinus reflex.
- ▶ The stimulation must not be applied on swollen, infected or inflamed areas or skin eruptions, such as phlebitis, thrombophlebitis, varicose veins, etc.
- ▶ The stimulation must not be applied on or in the vicinity of cancerous cells.
- ▶ The stimulation must not be used in cases of numbness/ sensory disturbances in the area of application.
- ▶ Do not use the device during pregnancy unless it is prescribed by your doctor.
- ▶ Do not use in cases of severe local inflammation.
- ▶ Do not use if you are using an intrauterine device (contraceptive coil) containing metal.
- ▶ Do not use in cases of higher-grade hemorrhoidal disease (when using a rectal probe).
- ▶ The device should not be used in cases of bladder overflow incontinence.
- ▶ In cases of fistulae or prolapsed uterus/vagina, the device must not be used.
- ▶ The device must not be used if there are infections in the vaginal or rectal area.
- ▶ The electrodes and probes should only be used by one person to avoid the risk of cross-infections.

- ▶ The device must not be used in the case of complete denervation of the pelvic floor.

SIDE EFFECTS

The side effects of the TENS/EMS units are generally mild, even with prolonged use. If disturbing effects occur, interrupt the treatment and have a physician evaluate the further use of the device.

Skin irritation and minor burns are possible side effects. If these occur, discontinue use and consult your doctor.

If muscle stimulation is too intense and prolonged, muscle pain in the sense of muscle soreness may occur. To avoid an increase in pain, shorter and less intense stimulation should be used, especially during the first treatment sessions. During intense higher-frequency stimulation (from approx. 15 Hz), persistent muscle tension may occur. If the electrodes are applied in the chest area, this can lead to disruption of respiratory activity during stimulation.

WARNINGS & SAFETY INSTRUCTIONS

The stimulation emitted by the device may be sufficient to cause an electric shock. Electric current of this magnitude must not flow through the rib cage (thorax) or through the chest and the upper back, as it can cause cardiac arrhythmia.

Do not place electrodes over the heart!

- ▶ Do not place electrodes on the front of the throat as spasm of the laryngeal and pharyngeal muscles may occur. Stimulation over the cervical cavity (throat area) may close the airway, make breathing difficult, and have negative effects on the heart rhythm or blood pressure.
- ▶ Do not place electrodes on your head or in places where the electric current can flow transcerebrally (through the head).
- ▶ If the therapy is ineffective or becomes uncomfortable, stimulation should be discontinued until its use is re-evaluated by a physician.
- ▶ The device does not have AP/APG protection. Do not use it in an explosive atmosphere or close to flammable mixtures.

- ▶ The device may only be used according to medical advice.
- ▶ Electrodes must not be placed over the eyes, the mouth or inside the body.
- ▶ The device is to be used only when awake.
- ▶ Do not use the device near other electrical equipment or power sources.
- ▶ If you experience pain, skin irritation, dizziness, difficulty breathing or other unusual symptoms while using the device, discontinue use and consult a physician immediately.
- ▶ Use only original accessories from the manufacturer with your device (see scope of delivery).
- ▶ Do not use this device for undiagnosed pain symptoms before consulting a physician.
- ▶ Do not use this device while driving, operating machinery, near water, or during activities where involuntary muscle contractions may expose the user to excessive risk of injury.
- ▶ Always turn off the unit before applying to or removing accessories from the skin or mucosa.
- ▶ Keep this device out of the reach of children.
- ▶ The device has no curative effect.
- ▶ Electronic monitoring devices (such as ECG monitors and ECG alarms) may not function properly when electrical stimulation is used.
- ▶ The device should be cleaned with soapy water on a clean cloth after each use.
- ▶ The device is an electronic product and must therefore not be disposed of with residual waste. Further information at: <https://axion.shop/en/pages/environment-and-disposal>
- ▶ For first-time users, muscle stimulation can be an unfamiliar sensation. It is recommended to start at a low intensity in a sitting position to familiarize yourself with the feeling before increasing the intensity.
- ▶ Consult a physician before using the device after giving birth.
- ▶ If an IUD is employed for contraception, consult a physician before use of the device.

- ▶ Please wait until your menstrual period is over, because during this time heavy use of the abdominal muscles is not recommended.
- ▶ If the battery fluid leaks and has come into contact with the skin, rinse the affected area thoroughly with plenty of clean water.
- ▶ If the battery fluid comes in contact with your eyes, avoid rubbing your eyes. Flush your eyes thoroughly with plenty of water and seek medical attention immediately.
- ▶ If you notice unusual heat, odors, or flames coming from the battery, stop using the device.
- ▶ Simultaneous connection of a patient to a high frequency surgical equipment and the stimulator may result in burns at the site of the stimulator electrodes and possible damage to the stimulator.
- ▶ Do not operate the product within 1 m of shortwave or microwave therapy equipment as it may produce instability in the stimulator output.

MAINTENANCE, TRANSPORT, STORAGE

- ▶ The device must not be exposed to large quantities of liquids/ water.
- ▶ After use, place the device back into the foam insert of the storage case for protection.
- ▶ Store the device in the storage case in a dry and cool place if you are not going to use it for a longer period of time.
- ▶ The fully charged battery has an operating time of 2 to 3 hours.
- ▶ If you find that the battery cannot be charged effectively, or the fully charged battery can only be used for a very short time, contact our service department.
- ▶ When not in use, store the device only with a fully charged battery. The device should be fully recharged after 6 months at the latest.

SCOPE OF DELIVERY

Your STIM-PRO X9B is delivered with the following components. All parts can be reordered separately from our online shop at www.axion.shop if needed.


Quantity	Product name	Order no.
1 x	TENS/EMS device	1090
1 x	USB charger	20075
2 x	4x Self-adhesive electrodes 50x50mm	20007
1 x	4x Self-adhesive electrodes 100x50mm	20021
2 x	2x Connection cable	20032
1 x	Storage case	109009
1 x	Instructions for Use	109001

EN

CHARGING THE DEVICE

Please charge the battery completely before using the device.

- ▶ Connect the small end of the USB cable to the device and the larger end to the power adapter.
- ▶ Plug the adapter into a standard power outlet.
- ▶ The device cannot be used during the charging process. The battery symbol flashes on the LCD display while the device is charging. Once charging is complete, a battery symbol with a fully filled bar will be displayed continuously.
- ▶ Remove the power plug from the socket and disconnect the USB cable from the device.

PREPARATION FOR TREATMENT

Skin preparation

The correct preparation of the skin before placement of the electrodes will prolong the life of the electrodes and reduce the risk of skin irritation. After you have connected the electrodes to the stimulator, proceed as follows:


- ▶ Determine the locations where the electrodes are to be placed.
- ▶ Clean the area with water and mild soap which does not contain any oils; do not use alcohol. Rinse the skin and dry it thoroughly.
- ▶ Trim excess body hair with scissors; do not shave.
- ▶ When removing the electrodes, always pull them off in the direction of hair growth.

Connecting and placing the electrodes

- ▶ Connect the plug connector to the electrodes before placing the electrodes on the skin.
- ▶ Plug the other end of the cable into the output socket of the device.
- ▶ Do not pull on the cables, but only on the plugs attached to the ends of the cables.
- ▶ Please position the electrodes on the relevant areas of the body as suggested in the brochure 'Suggestions for electrode placement'.
Note: We remind you that it is recommended to consult a health professional to test the positioning of the electrodes for the best possible pain relief.
- ▶ Select a TENS or EMS mode.
- ▶ You may need help attaching the electrodes to areas that are hard to reach (lower and upper back).
- ▶ When stimulating the muscles of the arms or legs, remember that muscle contraction can cause involuntary limb movements that could injure you or others.
- ▶ Do not switch on the device until all electrodes or probes and cables are properly attached.

OPERATING THE DEVICE


- CH1 Setting Channel 1
- CH2 Setting Channel 2
- CH3 Setting Channel 3
- CH4 Setting Channel 4
- ⋮
|| All channels/pause function:
 - ▶ Press and hold for simultaneous setting of all 4 channels
 - ▶ Press briefly to pause
- OFF
- OFF and minus key:
 - ▶ Press and hold to switch off device
 - ▶ Press to decrease the parameters
- SET Parameter selection
- ON
+ ON and Plus key:
 - ▶ Press and hold to switch on the device
 - ▶ Press to increase the parameters

EN

1. Charging the device



- ▶ Charge the device fully before use.

EN


2. Connecting the electrodes/ probes

- ▶ Connect the device cables to the electrodes or probes and place them (see page 42 for electrodes, page 55 for probes).

3. Switching on the device

- ▶ Press and hold the ON+ key for two (2) seconds to switch on the unit.


4. Selecting preset programs

The device offers the option to select different programs for CH1-2 and CH3-4 or one program for all channels. Proceed as follows:

- ▶ To select a program for CH1 and CH2 only, press the SET key. The program display starts to flash. Press CH1 or CH2 and use ON/+ or OFF/- to select the desired program. The therapy type TENS, EMS or MASSAGE will be displayed depending on the selected program.
- ▶ To select a program for CH3 and CH4 only, press the SET key. The program display starts to flash. Press CH3 or CH4 and the intensity for CH3 and CH4 will be displayed. Press ON/+ or OFF/- to select the desired program. The therapy type TENS, EMS or MASSAGE will be displayed depending on the selected program.
- ▶ To select a program for all channels (CH1-CH4), press the SET key. The program display starts to flash. Select the desired program with the ON/+ or OFF/- key. The therapy type TENS, EMS or MASSAGE will be displayed depending on the selected program.

You can find the program overview on pages 50-58.



5. Setting the treatment time

- ▶ Once the program is selected, press the SET key again. The program stops flashing and the timer starts to flash. Press ON/+ or OFF/- to select the treatment duration (if available). The timer can be set from 5 minutes to 60 minutes or to continuous operation C in 5-minute increments. Once the desired treatment duration is selected, press SET again to return to the main screen.

EN


Note: When using different programs for CH1/2 and CH3/4, the treatment times can be set differently.

6. Setting the intensity/current


- ▶ The intensity can be adjusted from 0–99 mA for each channel. Select the channel you want to set by pressing the CH1, CH2, CH3 or CH4 key. CH1, CH2, CH3 or CH4 flashes on the display.
- ▶ To increase or decrease the intensity, press ON + (to increase) or OFF – (to decrease) repeatedly until the desired intensity level flashes on the display. Press and hold the ON/+ or OFF/– key to continuously increase or decrease the value.
- ▶ To adjust the intensity of all four channels at once, press the key  for two (2) seconds; all channels CH1, CH2, CH3, CH4 will flash on the display. Press ON + (to increase) or OFF – (to decrease) repeatedly until the desired intensity level flashes on the display.
- ▶ Press SET to save your selection.

EN

7. Switching off the device

- ▶ The device switches off automatically after the therapy session time has expired.
- ▶ To turn the device off manually, press and hold the OFF key for three (3) seconds. The display fades out and the device switches off.
- ▶ In an emergency, you can also disconnect the plug(s) from the device and then remove the electrodes.
- ▶ To prevent unpleasant electric shocks, never remove the electrodes while the device is still switched on.

ADDITIONAL FUNCTIONS

Pause function

At any time during stimulation, press the Pause key  to pause the impulse output and press it again to resume stimulation.

Key lock

Press and hold the ON + and OFF – keys simultaneously for 1 second to lock or unlock the device; the lock icon appears on the display.

Note: The pause function can still be used at any time if one of the intensities is greater than or equal to 1, even when the key lock is active.

Treatment time

- ▶ The timer can be set from 5 minutes to 60 minutes or to continuous operation C in 5-minute increments. The time is counted down on the display in 1-minute decrements for the duration of the session.
- ▶ The device switches off automatically once the therapy session time has expired.
- ▶ The last therapy time set is saved.

- ▶ If you change the program mode during therapy, the therapy time is not restarted unless you reset it.

Resetting the intensity level

For your safety and comfort, the intensity level is reset to 0 every time the device is turned off and at the end of a therapy session.

The treatment will be stopped if the electrodes are not placed correctly. If you change the mode during the therapy session, the device switches to the lowest intensity level and displays 0 on the screen.

EN

Battery indicator

The battery indicator displays low battery level indication when the remaining battery capacity reaches 20%. Please recharge the battery soon.

Treatment time memory

- ▶ To access the treatment time memory, turn on the device and press and hold the ON/+ and SET keys for 3 seconds. The display now shows the total application time.
- ▶ The last 60 sessions are recorded. In the treatment time memory, press the ON/+ or OFF/- button to navigate through the session numbers for each mode. Use the SET button to switch between the TENS, EMS, and Massage modes. The duration of the respective session will then be displayed. Press the CH1, CH2, CH3, or CH4 buttons to switch between channels.
- ▶ To delete a specific treatment time press and hold the SET key for 3 seconds.
- ▶ To delete all recordings press and hold the SET and OFF/- keys for 3 seconds.
- ▶ Press and hold the ON/+ and SET keys for 3 seconds to exit the recording mode.

TENS STIMULATION

Programs 1–16 have preset parameters which have proven beneficial for their respective fields of application. These settings may also be used for other indications under certain circumstances, as each patient responds in a highly individual way to stimulation.

EN

Always start with the lowest intensity and increase it until a tingling sensation is felt. All programs are different and therefore feel different. You can try all the programs at the start and choose the one that feels comfortable. Increase the intensity slowly and never to the point where it hurts. The stimulation should not feel uncomfortable. Start with short sessions of 5 or 10 minutes until your body gets accustomed to the stimulation.

Note: For all TENS programs, the standard application duration is 30 to 45 minutes.

TENS-programs P 01–P 16

P 01

Classic TENS mode

Pulse frequency: 80 Hz | Pulse width: 180 μ s

Applicable to most acute and chronic types of pain.
Gate control effect

P 02

Low frequency burst TENS

Pulse frequency: 2–100 Hz | Pulse width: 180 μ s

Endorphin release. For application on arms, legs, knees and lower back

P 03

Modulated TENS

Pulse frequency: 80 Hz | Pulse width: 70–180 μ s

To counteract habituation effect with prolonged use.
Applicable to most acute and chronic types of pain

P 04

Alternating frequency 15/2 Hz

Pulse frequency: 2–15 Hz | Pulse width: 180 μ s

Application on joints, neck and lower back

P 05

Alternating frequency 80/2 Hz

Pulse frequency: 2–80 Hz | Pulse width: 180 μ s

For long-term treatments. Effective for most pain indications

- P 06** **Trigger Point TENS**
Pulse frequency: 10 Hz | Pulse width: 180 μ s
Application for seasickness/travel sickness, nausea, and vomiting by placing the electrodes from behind around the cervical vertebra C6.
- P 07** **Migraine and sensitive areas of the body**
Pulse frequency: 80 Hz | Pulse width: 60 μ s
Due to the reduced pulse duration, ideal for applications on sensitive areas of the body, e.g., on the head or face
- P 08** **Chronic pain**
Pulse frequency: 3-10 Hz | Pulse width: 200 μ s
Use e.g. for pain in the neck, hips, shoulders, elbows, lumbago, menstrual cramps, knee osteoarthritis
- P 09** **Rheumatic pain**
Pulse frequency: 5-50 Hz | Pulse width: 25-250 μ s
Application on arms, legs, feet, also suitable for prolonged use thanks to modulation
- P 10** **Sweep mode**
Pulse frequency: 5-125 Hz | Pulse width: 120 μ s
In case of resistance to therapy, sensory and motor stimulation, anti-habituation
- P 11** **Frequency change & pulse duration modulation**
Pulse frequency: 7-100 Hz | Pulse width: 237-260 μ s
Can be used for most acute and chronic types of pain
- P 12** **Frequency modulation**
Pulse frequency: 7-80 Hz | Pulse width: 260 μ s
Use for pain in the back, knee, shoulder, leg, elbow and neck; rheumatic pain
- P 13** **Prof. Han stimulation**
Pulse frequency: 2-100 Hz | Pulse width: 200-300 μ s
Muscular spasms, allodynia. Duration 30 minutes
- P 14** **High frequency TENS**
Pulse frequency: 100 Hz | Pulse width: 250 μ s
Acute pain. Tension. Gate control effect

P 15 **Combination TENS mode**
 Pulse frequency: 3-100 Hz | Pulse width: 180-260 μ s

4-phase TENS. For intractable and acute pain. Includes high and low frequency TENS and modulation of frequency, pulse duration, and intensity. Duration 40 minutes

P 16 **Frequency sweep mode**
 Pulse frequency: 5-150 Hz | Pulse width: 120 μ s

Application on wrist and fingers, long-term use. Duration 40 minutes.

EMS STIMULATION

Please note that the following programs consist of several phases. Programs 17–28 contain a warm-up phase of 60 seconds at the beginning. This is followed by the actual stimulation, which can consist of different phases. The goal is to contract a muscle and then relax it. All programs cause contraction and differ mainly in the speed and duration of contractions. As with any exercise program, start slowly at a low intensity (5 to 10 minutes) to warm up. Depending on the progress of your muscle performance, you can increase the intensity and duration of treatment.

Use the device regularly for a long period of time to maintain the benefits you have achieved during the workout.

After strenuous exercise or physical exertion, always use a lower intensity to prevent muscle fatigue.

Note: For all EMS programs, the standard application duration is 5 to 20 minutes.

EMS-Programs P 17–P 33

P 17 **EMS – Strength (Beginner)**
 Pulse frequency: 10-40 Hz | Pulse width: 150-300 μ s

For small muscle groups

P 18 **EMS – Strength (Beginner)**
 Pulse frequency: 10-50 Hz | Pulse width: 150-350 μ s

For large muscle groups

P 19 EMS – Strength (Advanced)
Pulse frequency: 10-70 Hz | Pulse width: 150-300 μ s

For small muscle groups

P 20 EMS – Strength (Advanced)
Pulse frequency: 10-70 Hz | Pulse width: 150-400 μ s

For large muscle groups

P 21 EMS – Muscle mass (Beginner)
Pulse frequency: 10-50 Hz | Pulse width: 150-300 μ s

For small muscle groups

P 22 EMS – Muscle mass (Beginner)
Pulse frequency: 10-50 Hz | Pulse width: 150-350 μ s

For large muscle groups

P 23 EMS – Muscle mass (Advanced)
Pulse frequency: 10-65 Hz | Pulse width: 150-350 μ s

For small muscle groups

P 24 EMS – Muscle mass (Advanced)
Pulse frequency: 10-65 Hz | Pulse width: 200-400 μ s

For large muscle groups

P 25 EMS – Stamina (Beginner)
Pulse frequency: 6-18 Hz | Pulse width: 150-300 μ s

For small muscle groups

P 26 EMS – Stamina (Beginner)
Pulse frequency: 6-18 Hz | Pulse width: 150-350 μ s

For large muscle groups

P 27 EMS – Stamina (Advanced)
Pulse frequency: 6-15 Hz | Pulse width: 150-350 μ s

For small muscle groups

P 28 EMS – Stamina (Advanced)
Pulse frequency: 6-15 Hz | Pulse width: 150-400 μ s

For large muscle groups

P 29 EMS – Relax I
Pulse frequency: 100 Hz | Pulse width: 150-200 μ s

Muscle regeneration, muscle relaxation – moderate intensity

EN	P 30	EMS – Relax II Pulse frequency: 80-100 Hz Pulse width: 200-300 μ s
		Muscle regeneration, muscle relaxation – moderate intensity
	P 31	EMS – Regeneration I Pulse frequency: 35 Hz Pulse width: 150 μ s
		For small muscle groups
	P 32	EMS – Regeneration II Pulse frequency: 35 Hz Pulse width: 300 μ s
		For large muscle groups
	P 33	EMS – Dynamic stimulation Pulse frequency: 35 Hz Pulse width: 250 μ s
		Dynamic stimulation of large and small muscles

What is considered a large or a small muscle?

Muscle groups are divided into large or small muscle strands. Below you can find some examples:

Large muscles:

- ▶ Thigh muscles
- ▶ Gluteal muscles
- ▶ Abdominal muscles
- ▶ Back muscles
- ▶ Chest muscles

Small muscles:

- ▶ Calf muscles
- ▶ Triceps, biceps
- ▶ Shoulder muscles
- ▶ Forearms
- ▶ Fingers, toes

Pelvic floor-programs P34 - P41

Note:

- ▶ These programs are for treatment with self-adhesive electrodes or vaginal/rectal electrodes.
- ▶ These programs involve noticeable but still comfortable contractions of the pelvic floor muscles.
- ▶ To ensure optimal stimulation, we recommend the use of contact gel with each application.
- ▶ The standard pelvic floor training session lasts 15 minutes and should be performed once per day.

Proceed as follows:

► **Plug in connecting cable**

Remove the probe from the packaging and clean it under running water. Connect the connection cable to the device and to the vaginal/rectal probe.

► **Insert the probe**

If necessary, visit the toilet before starting the procedure. The procedure should be done lying down with slightly drawn knees. Apply contact gel to the probe before insertion. Insert the probe slowly and gently into the vagina or anus.

► **Select the program and start the stimulation**

Switch on the device. Select the desired program (P34–P41) and start the stimulation by increasing the current intensity with the ON/+ and OFF/- keys of the relevant channel.

The program starts immediately after the intensity is set!

► **End the stimulation**

Switch off the device. Pull the probe out of the vagina or anus slowly.

Caution: To avoid breaking the cable, the probe must always be pulled out by the shaft! Clean the probe with lukewarm soapy water and rinse it.

Note: Please use the programs in daily rotation

P 34	Stress-incontinence	Pulse frequency: 50 Hz Pulse width: 300 μ s
P 35		Pulse frequency: 45 Hz Pulse width: 200 μ s
P 36		Pulse frequency: 45-50 Hz Pulse width: 200-300 μ s

A common cause of stress incontinence is a defect of the closure apparatus, which usually results from untrained or overstretched pelvic floor muscles, for example after heavy childbirth. Age-related hormone deficiency in women may also play a role.

P 37 **Strengthening of the pelvic floor**
Pulse frequency: 35 Hz | Pulse width: 250 μ s

For maintenance or general training of the pelvic floor muscles. For muscle maintenance of the pelvic floor, one workout per week is sufficient.

Note: Please use the programs in daily rotation

P 38	Urge-	Pulse frequency: 10 Hz Pulse width: 200 μ s
-------------	--------------	---

P 39	incontinence	Pulse frequency: 15 Hz Pulse width: 200 μ s
-------------	---------------------	---

With urge incontinence, there is often an intense urge to urinate, even with small amounts of urine, which can lead to premature urination. This type of incontinence is caused by overactive bladder muscles.

Note: Please use the programs in daily rotation

P 40	Mixed-	Pulse frequency: 10-50 Hz Pulse width: 200-300 μ s
-------------	---------------	--

P 41	Incontinence	Pulse frequency: 10-50 Hz Pulse width: 200-300 μ s
-------------	---------------------	--

Mixed incontinence is a form of urinary incontinence in which symptoms of both stress incontinence and urge incontinence occur. This means that both involuntary leakage of urine during physical exertion or pressure on the bladder and a strong urge to urinate accompanied by sudden leakage of urine may occur despite a small amount of urine.

EMS Massage programs P 42–P 57

An EMS massage uses electrical impulses to cause slight muscle contractions. This type of EMS application is intended to contribute to muscle relaxation (**scientifically investigated in several clinical studies, e.g. by Hsueh TC, Cheng PT, Kuan TS, Hong CZ**).^{*} The intensity is slowly increased until the muscle begins to contract

Note: For all Massage programs, the standard application duration is 15 to 30 minutes.

P 42	Massage	Pulse frequency: 1,25-90 Hz Pulse width: 90 μ s
-------------	----------------	---

Lower back

P 43	Massage	Pulse frequency: 3,3-58 Hz Pulse width: 90 μ s
-------------	----------------	--

Lower back, deltoid, elbow, knee, abdominal muscles, shoulder, sciatica, hip and thigh, calf, ankle/foot

P 44	Massage	Pulse frequency: 1,25-90 Hz Pulse width: 90 μ s
-------------	----------------	---

Lower back, deltoid, elbow, knee, abdominal muscles, shoulder, sciatica, upper arm, hip and thigh, calf, ankle/foot

- P 45** **Massage** Pulse frequency: 12,5-90 Hz | Pulse width: 90 μ s
Lower back, deltoid, knee, abdominal muscles, shoulder, sciatica, hip and thigh, ankle/foot
- P 46** **Massage** Pulse frequency: 3,3-90 Hz | Pulse width: 90 μ s
Lower back, deltoid, elbow, knee, sciatica, upper arm, hip and thigh, calf, ankle/foot
- P 47** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Lower back, deltoid, knee, sciatica, upper arm, hip and thigh, calf
- P 48** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Lower back, knee, sciatica, calf
- P 49** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Deltoid
- P 50** **Massage** Pulse frequency: 3,30-58 Hz | Pulse width: 90 μ s
Knee
- P 51** **Massage** Pulse frequency: 3,30-90 Hz | Pulse width: 90 μ s
Lower back
- P 52** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Lower back, deltoid, abdominal muscles, shoulder, sciatica, hip and thigh, ankle/foot
- P 53** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Elbow, knee, abdominal muscles, shoulder, calf
- P 54** **Massage** Pulse frequency: 3,30-58 Hz | Pulse width: 90 μ s
Elbow, knee, abdominal muscles, shoulder, calf
- P 55** **Massage** Pulse frequency: 3,30-90 Hz | Pulse width: 90 μ s
Deltoid, elbow, ankle/foot
- P 56** **Massage** Pulse frequency: 1,60 Hz | Pulse width: 90 μ s
Deltoid
- P 57** **Massage** Pulse frequency: 1,25-90 Hz | Pulse width: 90 μ s
Deltoid, elbow, abdominal muscles, shoulder, sciatica

*Hsueh TC, Cheng PT, Kuan TS, Hong CZ. The immediate effectiveness of electrical nerve stimulation and electrical muscle stimulation on myofascial trigger points. Am J Phys Med Rehabil. 1997;76(6):471-6.

MANUAL MODE

Manual Programs P 58–P 60

P 58	EMS Mode
-------------	-----------------

freely adjustable parameters

P 59	TENS Mode
-------------	------------------

freely adjustable parameters

P 60	TENS Burst Mode
-------------	------------------------

freely adjustable parameters (fixed frequency 100 Hz)

The device offers the option to use free parameters in TENS or EMS mode. To do this, proceed as follows:

TENS

1. Therapy mode selection (TENS)

Follow step 4 on page 45 to select program P 59 or 60 for manual TENS mode.

2. Setting the pulse width (Width)

Press the SET key; WIDTH starts blinking; set the pulse width with the ON/+ or OFF/- key. The pulse width can be adjusted from 10 to 400 μ s. If there are no specific instructions from the therapist, select the range from 70 to 120 μ s.

3. Setting the frequency (Rate)

Press SET again; RATE (frequency) starts blinking on the display; set the frequency with the ON/+ or OFF/- key. The frequency is adjustable from 1 Hz to 150 Hz. If the therapist has not specified the frequency for the stimulation, please use the range between 70 Hz and 120 Hz

4. Setting the treatment time

Press SET again; the treatment time setting flashes; set the treatment time with the ON/+ or OFF/- key.

Note: For all TENS treatments, the standard application duration is 30 to 45 minutes.

5. Setting the intensity/current

Follow point 6 "Setting the intensity/current" on page 47 and increase it with the ON/+ or OFF/- key.

EMS

1. Therapy mode selection (EMS)

Follow step 4 on page 45 to select program P 58 for manual EMS mode.

2. Setting the pulse width (Width)

Press the SET key; WIDTH starts blinking; set the pulse width with the ON/+ or OFF/- key. The pulse width can be adjusted from 10 μ s to 400 μ s. If there are no specific instructions from the therapist, select the range from 70 μ s to 120 μ s.

3. Setting the frequency (Rate)

Press SET again; RATE (frequency) starts blinking on the display; set the frequency with the ON/+ or OFF/- key. The frequency is adjustable from 1 Hz to 150 Hz. If the therapist has not specified the frequency for the stimulation, please use the range between 70 Hz and 120 Hz.

4. Setting the Ramp

Press SET again; RAMP starts blinking on the display. Use the ON/+ or OFF/- key to set the ramp time from 1–8 seconds. The term "Ramp" refers to the function used to increase (ramp up) or decrease (ramp down) the intensity of the electric current to gently activate the muscle.

- ▶ Press SET again; RAMP UP starts blinking; set the time for the increased ramp time with the ON/+ or OFF/- key.
- ▶ Press SET again; RAMP DOWN starts blinking on the display; set the time for the decreased ramp time with the ON/+ or OFF/- key.

5. Setting the contraction time (ON time)

Press SET again; ON time starts blinking, which determines the duration of the contraction time. By pressing the ON/+ or OFF/- key, the contraction time can be set from 0–60 seconds. The total signal cycle includes both ramp times (increase and decrease).

Example: 3 seconds ramp-up time + 10 seconds ON time + 4 seconds ramp-down time = 17 seconds total cycle time

6. Setting pause times (OFF time)

Press SET again; OFF time starts blinking, which determines the duration of the pause between cycles. By pressing the ON/+ or OFF/- key, the pause time can be set from 0–60 seconds.

7. Setting the treatment time


Press SET again; the treatment time setting flashes; set the treatment time with the ON/+ or OFF/- key

Note: For all EMS treatments, the standard application duration is 5 to 20 minutes.

8. Setting the intensity/current

Follow point 6 “Setting the intensity/current” on page 47 and increase it with the ON/+ or OFF/- key.

LABEL

EN

	Manufacturer	REF	Article number
CE 2460	Conformity marking with number of the notified body (device complies with the requirements of (EU)2017/745)	MD	Medical device
		UDI	Unique Device Identification
	Date of manufacture YYYY-MM-DD	SN	Serial number
	Caution! Check the instructions for use for important safety-related information such as warnings and precautions.		Follow instructions for use
			Follow operating instructions
IP22	IP protection class		Temperature limit
	Device type BF		Temperature limit for storage/transportation
	Protect from sunlight		Protect from moisture
	Do not dispose of with residual waste		2D datamatrix barcode: encoded EAN, date of manufacture and serial number

FUNCTION AND SAFETY TEST

For safety reasons, the following precautions should be carried out once a week.

- ▶ Check the device for external damage: deformation of the housing or damage to the electrode output sockets.
- ▶ Check that no type plates or stickers are damaged.
- ▶ Check cables and electrodes for damages.
- ▶ A safety inspection and maintenance must be carried out by authorized personnel before the device is used for the first time and before operating it after a period of inactivity; at a minimum every 24 months.
- ▶ The instructions for use must always accompany the unit. If there are any problems or defects, please contact your retailer.

REPAIRING MINOR FAULTS

If the device does not work properly, the instructions below may help identify possible faults. If the problem cannot be solved with any of these tips, contact the nearest service center.

- ▶ Check cables and connections. Have all plugs been connected correctly?
- ▶ If a cable is damaged externally, replace it immediately with an undamaged one. Only use original accessories.
- ▶ No current is flowing? Check if the electrodes are correctly connected to the same channel and if the current is strong enough.

TECHNICAL DESCRIPTION

01	Channels	4 channels, each with adjustable intensity
02	Pulse form	Symmetrical, bi-phasic square pulse
03	Product specification	TENS/EMS/MASSAGE
04	Intensity	Adjustable, 0–99 mA with a load of a minimum of 500 Ohm per channel
05	USB charging port	USB Type C
06	Voltage	50 V
07	Power source	Rechargeable lithium-ion battery
08	Battery indicator	Low battery level indication at 20% remaining capacity.
09	Timer	5–60 minutes (in 5-minute increments) and continuous operation C
10	Display	LCD with backlighting
11	Size	13.7 cm (L) x 7.8 cm (W) x 2.3 cm (H)
12	Weight	170 g incl. battery
13	Frequency	Adjustable from 1–150 Hz
14	Active time (ON time)	Adjustable from 1–60 seconds
15	Pause time (OFF time)	Adjustable from 0–60 seconds
16	Pulse width	Adjustable from 10–400 μ s
17	Ramp	Adjustable from 1–8 seconds, in 1-second increments
18	Function mode	TENS: 18 programs EMS: 26 programs MASSAGE: 16 programs
19	Operating conditions	Temperature: 10° C to 40° C Relative humidity: 40–90% Air pressure: 700–1013 hPa Altitude in metres: 3000
20	Transport and storage conditions	Temperature: -10° C to 60° C Relative humidity: 30–95% Air pressure: 500–1060 hPa
21	Notice	All technical values include a tolerance of +/- 5%

WARRANTY

All devices come with a warranty of 24 months from the date of purchase. The warranty applies to the device only (except for the battery, which is covered by a limited warranty of 6 months) and includes spare parts and labor costs. The warranty does not apply to damage resulting from improper handling, the failure to follow the instructions for use, loss or dropping.

EN

CONFORMITY TO SAFETY STANDARDS

The STIM-PRO X9B device is in compliance with the safety standards EN 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance & EN 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances - Requirements and tests.

Special precautions, information and characteristics on “Electromagnetic compatibility (EMC)”: Medical electrical devices are subject to special requirements and precautions with regard to electromagnetic compatibility.

Warnings:

- ▶ The use of accessories, transducers and wiring other than those specified or provided by axion of this appliance may result in increased ELECTROMAGNETIC RADIATION or decreased electromagnetic immunity of the appliance and result in improper operation.
- ▶ Portable RF communications equipment (radios), including accessories such as antenna cables and external antennas, should be used no closer than 30 cm (12 inches) to any part or line of the TENS/EMS DEVICE as specified by the manufacturer. Failure to do so may result in a reduction in the performance characteristics of the device.
- ▶ Use of this appliance adjacent to other appliances or stacked with other appliances should be avoided as this may result in improper operation. If use in the manner described above is nevertheless necessary, this appliance and the other appliances should be observed to ensure that they are operating correctly.

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axion GmbH